

About us



0.9
ML
0.5


Maxxline is an affiliated brand of AHN Biotechnologie GmbH, which has acquired extensive experience and expertise in the production of plastic consumables for laboratories since 1999.

The market demands products that ensure and optimize the productivity and provide greater process reliability.

Maxxline reflect innovative “ways and means” for the production of plastic products, combined with a focus on the highest level of quality.

Based on this we have developed a new unique product - Maxxline Cryovial - with a StraightLine™ design.

Maxxline products are designed and manufactured in Germany.



Our star



The graphic symbol in our Maxxline logo is a Sea Star.


The Sea Star symbolizes our unique products.

At first sight, the Sea Star seems simple, and yet they are complex multi-cellular organisms with unique and surprising properties. The Sea Star, for example, can regenerate lost body parts.

Maxxline Cryovial is a new product, with surprising properties. A unique star when it comes to the storing of biological material, human or animal cells for research.

Like the Sea Star with its five arms, our Maxxline Cryovials stand solid on five feet.

Your tube, your future.




We make
the difference



As a leading manufacturer of cryogenic tubes we know that functionality and quality are properties of the highest importance in a product for professional use. Properties we always take into consideration when designing new products for our range.

It is important to us that our customer has the best experience when using our products. Therefore, convenience is a milestone in everything we do.

The development and final design of our Cryovials is based on many years of experience from working with and selling to leading laboratory technicians all over the world.



Maxx features



Ergonomics


Maxxline Cryovials have a single-handed quick seal closure. The super-fast thread design allows tightening or removal with only half a turn, avoiding long-term injuries to the user.

Max volume

Maxxline Cryovials are easy to read. The red line indicator shows the maximum recommended filling level - your security for safe storage.

Contamination

The StraightLine™ design of the Maxxline Cryovial gives you approximately 23% additional storage capacity compared to conventional cryo solutions.



Think ahead


We work exclusively with the best materials and product designs that optimize manufacturing.

The use of materials and power consumption are kept to a minimum.

Our products are produced of raw materials which are compliant to applicable regulations and standards, and are free of any substance considered dangerous for the environment.

Quality management system according to ISO 9001:2015.

Your tube, your choice, your possibilities, your future.





AHN Biotechnologie GmbH
Uthleber Weg 14 | 99734 Nordhausen | Germany
Phone: +49(0)3631/65242-0
E-mail: info@cappahn.com | Website: www.cappahn.com