

Plastic Cryogenic Boxes for 15mL / 50mL tubes

Maxxline centrifuge tube storage boxes are uniquely tailored to meet the demands of storing 15mL and 50mL conical tubes. These racks have a wide scope of applications, as described below.

Durable Polypropylene Storage Boxes

Maxxline centrifuge tube storage boxes are manufactured from a customized blend of extremely durable polypropylene durable. For this reason, they safely hold 15mL and 50mL conical tubes at sample freezing (-20°C) and cryopreservation (-90°C) temperatures, and can be sterilized at +121°C without losing structural integrity.

Vented Conical Tube Storage Boxes

To increase their thermal efficiency, Maxxline Cryogenic Boxes feature bottom-facing ventilation holes, which greatly increase the rate of freezing. These downward facing vents allow for minimal fluid retention.

Multiple Layout Options

Maxxline conical tube boxes accommodate up to 25 tubes in two different grid layouts. The 5x5 grid (CM25-125(X)) has a regular pattern that can hold up to twenty-five 15mL tubes. Larger conical tubes would fit in the 12-place layout that holds up to ten 50mL tubes and two 15mL tubes.

Multi-purpose Conical Tube Box

Maxxline centrifuge tube storage boxes can be adapted to various uses in the lab. This includes freezing 15mL and 50mL tubes as well as sorting and holding samples during centrifugation. Maxxline storage boxes may also be used for transporting samples in a dependable manner.

Highlights:

- A snap on lid closure secures samples in each box
- Vent holes allow for rapid freezing of samples stored on each rack and hassle-free liquid drainage
- Stock footprint (130x130mm) that fits well in most freezer racks
- Optional 12 or 25 place dividers hold many sample tubes on each rack
- Clear plastic construction improves sample visibility

Ordering Information

Cat. No.	Description	Dimensions	Colour	Packaging
CM12-125N	Plastic box, 130mm, grid for 10x50mL and 2x15mL tubes	130x130x130mm	Transparent	36 pcs.
CM12-125B	Plastic box, 130mm, grid for 10x50mL and 2x15mL tubes	130x130x130mm	Blue	36 pcs.
CM12-125R	Plastic box, 130mm, grid for 10x50mL and 2x15mL tubes	130x130x130mm	Red	36 pcs.
CM12-125G	Plastic box, 130mm, grid for 10x50mL and 2x15mL tubes	130x130x130mm	Green	36 pcs.
CM12-125Y	Plastic box, 130mm, grid for 10x50mL and 2x15mL tubes	130x130x130mm	Yellow	36 pcs.
CM25-125N	Plastic box, 130mm, grid for 25x15mL tubes	130x130x130mm	Transparent	36 pcs.
CM25-125B	Plastic box, 130mm, grid for 25x15mL tubes	130x130x130mm	Blue	36 pcs.
CM25-125R	Plastic box, 130mm, grid for 25x15mL tubes	130x130x130mm	Red	36 pcs.
CM25-125G	Plastic box, 130mm, grid for 25x15mL tubes	130x130x130mm	Green	36 pcs.
CM25-125Y	Plastic box, 130mm, grid for 25x15mL tubes	130x130x130mm	Yellow	36 pcs.
CM40N	Plastic box, 50mm, no dividers	130x130x50mm	Transparent	36 pcs.
CM40N-80	Plastic box, 80mm, no dividers	130x130x80mm	Transparent	36 pcs.
CM40	Plastic box bottom only, 50mm, no dividers	130x130x50mm	Transparent	36 pcs.